

SkinAxis, LLC CCIT, 675 US Highway One North Brunswick NJ, 08902-3378

Tel: 347-707-0829 Fax: 732-745-7270 www.skinaxis.com

Human Macrophages

Product Description

Human macrophages differentiated from Peripheral Blood Mononuclear Cells (PBMC)-derived monocytes.

Sample Collection and Processing


CD14⁺ monocytes were isolated from buffy coats of normal human volunteers by positive selection using an immunomagnetic cell separation method. Monocytes were differentiated into macrophages by culturing cells for 10-13 days in the SkinAxis proprietary medium.

Purity

The purity of PBMC-derived CD14⁺ monocytes, as measured by flow cytometric analysis, is ≥95% (left panel). Macrophages are >99% positive for the pan-macrophage CD68 marker and show efficient phagocytic function (right panel).

FACS analysis of CD14⁺ cells before (Total PBMC) and after purification indicates the enrichment of CD14⁺ monocytes population to >95%.

Purity and phagocytic function of macrophages



Staining of monocyte-derived macrophages with CD68 (left panel). The phagocytosis of fluorescent polystyrene beads confirms efficient phagocytic function of macrophages (right panel).

Thawing and General Culture Protocol

- 1. Wipe the frozen cryovial with 70% ethanol in the biosafety sterile hood. Partially thaw the vial of cells in the 37°C water bath for 2 mins. Thawing the cells at 37°C for longer than 2 minutes results in reduced viability. DO NOT submerge the entire vial in the water bath. Wipe the cryovial with 70% ethanol and move to a biosafety sterile hood.
- 2. Carefully add 500ul of DMEM supplemented with 10% FBS, 2mM Glutamine (complete medium) into the vial and quickly transfer the whole volume into the 10 ml of complete medium. DO NOT pipet repeatedly.
- 3. Plating at a density of 10⁵ cells per cm² is recommended. Incubate at 37°C, 5% CO₂ for 16-24 hours in a humidified incubator.
- 4. After 16-24 hours of culturing, check the cells under a light microscope; macrophages should have adhered to the culture dish/plate. Change the medium to remove any residual DMSO from cryopreservation.
- 5. Replace medium every 2 days.

Stability and Storage

Macrophages are cryopreserved in the SkinAxis Freezing Medium and are stable for long-term storage at -152°C or in liquid nitrogen. Storage at -80°C is NOT recommended. Once thawed, the cells must be used immediately.

Quality control

Cells are derived from donors with undetectable HIV-1, hepatitis B and hepatitis C. Macrophages are characterized by morphological analysis, specific marker expression, and phagocytic ability.

Handling

Primary cells from human blood are a potential biohazard. Treat as potentially infectious. All human sourced products should be handled at the Biological Safety Level 2 to minimize exposure of potentially infectious products.

THIS PRODUCT IS FOR IN VITRO RESEARCH USE ONLY